# Laboratory Preparation for 2019 Novel Coronavirus

National Influenza Centre of Myanmar (NIC)

National Health Laboratory

Date: 24.1.2020



# Overview

- Coronavirus
- Sample collection and shipment
- Lab Diagnosis
- Constraints

# Coronavirus

---

# Coronavirus

- Coronaviruses are named for the crown-like spikes on their surface
- Family *Coronaviridae*, subfamily *Coronavirinae*
- There are **four genera** - alpha, beta, gamma, and delta
- High frequency of recombination
- Common in a wide range of mammalian and avian species, most notably bats
- Human coronaviruses were first identified in the mid-1960s

# Seven coronaviruses that can infect people

## ➤ Common human coronaviruses

1. 229E (alpha coronavirus)
2. NL63 (alpha coronavirus)
3. OC43 (beta coronavirus)
4. HKU1 (beta coronavirus)

## ➤ Other human coronaviruses

5. MERS-CoV (the beta coronavirus that causes Middle East Respiratory Syndrome).
6. SARS-CoV (the beta coronavirus that causes severe acute respiratory syndrome).
7. 2019 Novel Coronavirus (2019-nCoV)

# 2019 Novel Coronavirus

- 2019- nCoV is a new strain that has not been previously identified in humans
- Zoonosis

# Sample collection and shipment

---

# Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases

## Interim guidance

### 14 January 2020

WHO/2019-nCoV/Laboratory/2020.2

## 1. Introduction

The purpose of this document is to provide interim guidance to laboratories and stakeholders involved in laboratory testing of patients who meet the definition of suspected case of pneumonia associated with a novel coronavirus identified in Wuhan, China (See: [Surveillance case definitions for human infection with novel coronavirus, interim guidance](#)).

Various existing WHO documents have been adapted for use in the drafting of this document, including WHO laboratory guidance for MERS-CoV (1-11). As information about the etiology, clinical manifestations and transmission of disease in the cluster of respiratory disease patients identified in Wuhan is evolving, WHO continues to monitor developments and will revise these recommendations as necessary.

The etiologic agent responsible for the cluster of pneumonia cases in Wuhan has been identified as a novel betacoronavirus, (in the same family as SARS-CoV and MERS-CoV) via next generation sequencing (NGS) from cultured virus or directly from samples received from several pneumonia patients. Electron microscopy revealed a virus with a characteristic crown morphology: a coronavirus. Working directly from sequence information, the team developed a series of genetic amplification (PCR) assays used by laboratories associated with the China CDC to detect several dozen cases as of today.

Full genome sequence data from the viruses have been shared officially with WHO and on the GISAID platform (<https://gisaid.org/>) and can inform the development of specific diagnostic tests for this emergent coronavirus. It is expected that validated PCR tests will become available soon. Until that time, the goals of diagnostic testing are to detect conventional causes of pneumonia early, to support disease control activities, and to work with reference laboratories that can perform pan coronavirus detection and directed sequencing.

## 2. Suspected case definition

For case definition see: [WHO Surveillance case definitions for human infection with novel coronavirus](#).

## 3. Specimen collection and shipment

Rapid collection and testing of appropriate specimens from suspected cases is a priority and should be guided by a laboratory expert. As extensive testing is still needed to confirm the 2019-nCoV and the role of mixed infection has not been verified, multiple tests may need to be performed and sampling sufficient clinical material is recommended. Local guidelines should be followed regarding patient or guardian's informed consent for specimen collection, testing and potentially future research.

Assure SOPs are available, and the appropriate staff is trained and available for appropriate collection, specimen storage, packaging and transport. There is still limited information on the risk posed by the reported coronavirus found in Wuhan, but it would appear samples prepared for molecular testing could be handled as would samples of suspected human influenza (2, 7-9). Attempts to culture the virus require a higher level of biosecurity.

**Samples to be collected (see Table 1 for details on sample collection and storage):**

1. Respiratory material\* (nasopharyngeal and oropharyngeal swab in ambulatory patients and sputum (if produced) and/or endotracheal aspirate or bronchoalveolar lavage in patients with more severe respiratory disease)
2. Serum for serological testing, acute sample and convalescent sample (this is additional to respiratory materials and can support the identification of the true agent, once serologic assay is available)

\*Modifiable with information on whether upper or lower respiratory material is better for coronavirus detection.

*A single negative test result, particularly if this is from an upper respiratory tract specimen, does not exclude infection. Repeat sampling and testing, lower respiratory specimen is strongly recommended in severe or progressive disease. A positive alternate pathogen does not necessarily rule out either, as little is yet known about the role of coinfections.*

Reference 2, 3, 7

Table 1. Specimens to be collected from symptomatic patients  
Guidance on specimen collection (adapted from reference 5)

| Specimen type   | Collection materials                                       | Transport to laboratory | Storage till testing                 | Comment  |
|---|--|-------------------------|--------------------------------------|--|
| Nasopharyngeal and oropharyngeal swab   | Dacron or polyester flocked swabs*                         | 4 °C                    | ≤5 days: 4 °C<br>>5 days: -70 °C     | The nasopharyngeal and oropharyngeal swabs should be placed in the same tube to increase the viral load. |
| Bronchoalveolar lavage  | sterile container *  | 4 °C                    | ≤48 hours: 4 °C<br>>48 hours: -70 °C | There may be some dilution of pathogen, but still a worthwhile specimen                                  |
| (Endo)tracheal aspirate, nasopharyngeal aspirate or nasal wash                | sterile container *  | 4 °C                    | ≤48 hours: 4 °C<br>>48 hours: -70 °C |  |
| Sputum  | sterile container  | 4 °C                    | ≤48 hours: 4 °C<br>>48 hours: -70 °C | Ensure the material is from the lower respiratory tract  |
| Tissue from biopsy or autopsy including from lung                             | sterile container with saline                              | 4 °C                    | ≤24 hours: 4 °C<br>>24 hours: -70 °C |  |
| Serum (2 samples acute and convalescent possibly 2-4 weeks after acute phase) | Serum separator tubes (adults: collect 3-5 ml whole blood) | 4 °C                    | ≤5 days: 4 °C<br>>5 days: -70 °C     | Collect paired samples:<br>• acute – first week of illness<br>• convalescent – 2 to 3 weeks later        |
| Whole blood   | collection tube  | 4 °C                    | ≤5 days: 4 °C<br>>5 days: -70 °C     | For antigen detection particularly in the first week of illness  |
| Urine   | urine collection container                                 | 4 °C                    | ≤5 days: 4 °C<br>>5 days: -70 °C     |  |

\*For transport of samples for viral detection, use VTM (viral transport medium) containing antifungal and antibiotic supplements. For bacterial or fungal culture, transport dry or in a very small amount of sterile water. Avoid repeated freezing and thawing of specimens.

Aside from specific collection materials indicated in the table also assure other materials and equipment are available: e.g. transport containers and specimen collection bags and packaging, coolers and cold packs or dry ice, sterile blood-drawing equipment (e.g. needles, syringes and tubes), labels and permanent markers, PPE, materials for decontamination of surfaces.

## Safety procedures during sample collection and transport

All specimens collected for laboratory investigations should be regarded as potentially infectious, and HCWs who collect, or transport clinical specimens should adhere rigorously to infection prevention and control guidelines and national or international regulations for the transport of dangerous goods (infectious substances) to minimize the possibility of exposure to pathogens (14). Implement the appropriate infection prevention and control precautions, guidance on IPC for the 2019-nCoV has been drafted (11).

## Assure good communication with the laboratory and provide needed information

To assure proper and fast processing of samples and to assure adequate biosafety measures in the laboratory, communication and information sharing is essential. Be sure you have alerted the laboratory of the urgency and situation before sending the sample. Also assure that specimens are

correctly labelled, and diagnostic request forms are filled out properly and clinical information is provided (see box: information to be recorded).

- Information to be recorded:**
- Patient information – name, date of birth, sex and residential address, unique identification number, other useful information (e.g. patient hospital number, surveillance identification number, name of hospital, hospital address, room number, physicians' name and contact information, name and address for report recipient).
  - Date and time of sample collection,
  - Anatomical site and location of specimen collection,
  - Tests requested.
  - Clinical symptoms and relevant patient history (including vaccination and antimicrobial therapies received, epidemiological information, risk factors).

# Samples to be collected

- Respiratory material (nasopharyngeal and oropharyngeal swab in ambulatory patients and sputum (if produced) and/or endotracheal aspirate or bronchoalveolar lavage in patients with more severe respiratory disease)
- Serum for serological testing, acute sample and convalescent sample (this is additional to respiratory materials and can support the identification of the true agent, once serologic assay is available)
  - *A single negative test result, particularly if this is from an upper respiratory tract specimen, does not exclude infection. Repeat sampling and testing, lower respiratory specimen is strongly recommended in severe or progressive disease. A positive alternate pathogen does not necessarily rule out either, as little is yet known about the role of coinfections.*

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- Already distributed to State and Regional Hospital Labs by email
- NHL FB page and Website

# 2019-nCoV Lab Request Form

**Novel Coronavirus (nCoV) Suspected Cases Laboratory request Form**

Please complete this form carefully and circle the response.

|  |                |   |                |
|--|----------------|---|----------------|
| <b>1. Report/Investigation Information:</b> Name of Investigator(s): _____                             |                |   |                |
| Name of Hospital: _____  |                |   |                |
| Date - Case Reported: ____/____/____   |                |   |                |
| Date - Case Investigated: ____/____/____   |                |   |                |
| <b>2. Case Identification:</b> Patient's Name: _____   |                |   |                |
| Date of Birth: ____/____/____ Age: years ____ months ____ Sex: ____                                    |                |   |                |
| Father's Name: _____ Mother's Name: _____  |                |   |                |
| Full Permanent Address: State/Region: _____ Township: _____  |                |   |                |
| Village/ward: _____ Street No. & House No: _____   |                |   |                |
| Phone No: _____  |                |   |                |
| <b>3. Hospitalization:</b> Yes / No  |                | Date of Hospitalization: ____/____/____ |                |
| Name of Hospital: _____  |                | Hospital Registration Number: _____     |                |
| Clinical Diagnosis: _____  |                |   |                |
| Outcome: Recovered completely / Death / Unknown  |                |   |                |
| <b>4. Travel History within last two weeks:</b> Yes / No   |                |   |                |
| Where: _____ When: ____/____/____  |                |   |                |
| <b>5. Poultry Contact History within last two weeks:</b> Yes / No                                      |                |   |                |
| <b>6. Sign and Symptoms:</b> Date of onset of first symptoms: ____/____/____                           |                |   |                |
| Fever: Yes / No / Unknown Cough: Yes / No / Unknown  |                |   |                |
| Lower respiratory tract involvement: dyspnea: Yes / No / Unknown or difficulty breathing: Yes / No     |                |   |                |
| Upper respiratory tract involvement: sore throat: Yes / No / Unknown or coryza: Yes / No               |                |   |                |
| Other symptoms: _____  |                |   |                |
| <b>7. Comorbid conditions:</b> Heart Disease ( ) Asthma ( ) Chronic Lung Disease ( ) Liver Disease ( ) |                |   |                |
| Immunocompromised ( ) Pregnancy ( ) Other (Specify): _____   |                |   |                |
| <b>8. Specimen Collection:</b>   |                |   |                |
|  | Date Collected | Date Sent to Lab                        | Date of Result |
| Nasopharyngeal Swab  | ____/____/____ | ____/____/____                          | ____/____/____ |
| Oropharyngeal Swab   | ____/____/____ | ____/____/____                          | ____/____/____ |
| Sputum   | ____/____/____ | ____/____/____                          | ____/____/____ |
| Bronchoalveolar Lavage   | ____/____/____ | ____/____/____                          | ____/____/____ |
| Endotracheal Aspirate  | ____/____/____ | ____/____/____                          | ____/____/____ |
| Serum  | ____/____/____ | ____/____/____                          | ____/____/____ |
| Whole blood  | ____/____/____ | ____/____/____                          | ____/____/____ |
| <b>9. Case Classification:</b> Lab confirmed Novel Coronavirus / Discard                               |                |   |                |
| <b>10. Signature of responsible person filling the form:</b> _____                                     |                |   |                |

- Already distributed to State and Regional Hospital Labs by email
- NHL FB page and Website

**STANDARD OPERATING PROCEDURES ON  
COLLECTION, STORAGE AND TRANSPORTATION FOR  
INFECTIOUS DISEASE SPECIMENS**

AND

**GUIDANCE ON LABORATORY SPECIMEN COLLECTION**

National Health Laboratory  
Department of Medical Services  
Ministry of Health and Sports

THE REPUBLIC OF THE UNION OF MYANMAR  
October, 2019 (Version 01)

Training Workshop on Infectious Substance Transport  
25th - 26th September 2019

# Receipt of Samples

- Assign NHL staff to provide 24-7 service for receipt of the samples
- Arrange the samples shipment to Reference lab as soon as possible

# Lab Diagnosis

---

- According to SEARO's suggestion, plan to send the suspected samples to National Institute of Health (NIH), Department of Medical Sciences, Thailand
- NIH, Thailand agrees to test suspected samples from Myanmar
- No charges to pay for testing
- WCO already discussed with World courier for shipment
- The shipment charges will be covered by WCO
- Turn around time – NIH said – results will be 2 days after the receipt of samples

# Constraints of NIC Myanmar

---

- Limited HR
- Require the particular PCR reagents and primers & probes for 2019- nCoV
- NIC is only ( enhanced ) Biosafety level 2 and handlings of novel coronavirus samples with current facilities are much be challenging to Biosafety and Biosecurity

## Expecting capacities in NIC, Myanmar

- Already have well-set molecular facilities in NIC, Myanmar
- Can perform in house lab diagnosis if specific reagents (primers and probes) with validated controls but require to send the samples to reference lab for confirmation
- Smart BSL 3 (2 rooms) are going to arise in NHL in 2022

# 23 January 2020 by WHO

- The Emergency Committee convened by the WHO Director-General under the International Health Regulations regarding the outbreak of novel coronavirus 2019 advised on 23 January 2020 that
  - The event did not constitute a Public Health Emergency of International Concern (PHEIC)
- Of 584 cases now reported, 575 of them and all 17 deaths have been in China

# References

- Diagnostic detection of Wuhan coronavirus 2019 by real-time RT-PCR, Berlin, 13.01.2020
- Interim guidance on Laboratory testing for 2019-nCoV in suspected human cases, 14.01.2020
- Coronavirus 101 presented by Respiratory Viruses Branch, Division of Viral Diseases, National Center for Immunizations and Respiratory Diseases, Centers for Disease Control and Prevention

**Thank you for your kind attention**